

Micro-Touch[®] NitraTex[®]



Powder Free Nitrile Examination Gloves
The synthetic choice with added protection

Ansell

Micro-Touch® NitraTex®

Powder Free Nitrile Examination Gloves

Micro-Touch NitraTex nitrile gloves provide a suitable alternative to natural rubber latex. Micro-Touch NitraTex offers high tear and chemical resistance. Available in two cuff lengths (240mm and 280mm) to offer extra protection in the wrist area.

Areas of use

- General wards.
- Scope cleaning (sleeve extenders available in polyethylene, vinyl or neoprene).
- Labs.
- Dental.
- Sterilising departments.

Glove colour

- Blue coloured nitrile glove for easy identification as a synthetic glove.

Special features

- Nitrile gloves offering a comfortable fit and good grip.
- Improved softer formula.
- Ambidextrous.
- Both powder and latex free.
- High tear resistance.
- Excellent resistance to a broad range of commonly used chemicals.
- Silicone free – ensures no delamination to contaminate test results.
- Thiurams and mercaptobenzothiazoles free – reduces the risk of contact dermatitis.
- Layer pack for easier dispensing and reduction of waste – save on glove wastage.
- Three year shelf-life.
- Size range: S-XL (four sizes).

Surface finish

- Mild textured finish for excellent grip.

Shape features

- Textured fingers and palm for excellent grip.

Cuff features

- Beaded cuff.
- Available in two cuff lengths: Standard Cuff (240mm) and Long Cuff (280mm) for extra protection when required.

Glove properties

- Material: 100 per cent nitrile.
- Colour: Blue.
- Donning agent: None.
- Design: Ambidextrous with beaded cuff.
- Fingertip thickness (single): 0.12mm.
- Cuff thickness (single): 0.08mm.
- Length: 240mm and 280mm.

Quality standards

- Statistically tested to an AQL of 1.5 for pinholes (exceeds world standards).
- Manufactured in accordance with ISO 13485.
- Carries the CE mark and conforms to EN 455 – 1, 2 & 3.
- FDA certified – Good Manufacture Practice (GMP).

Ordering information

Product	Code	Size	Packaging
Micro-Touch NitraTex (240mm cuff length)	4470	Small	100 gloves per box/ 10 boxes per carton
	4472	Medium	
	4474	Large	
	4476	Extra Large	
Micro-Touch NitraTex EP (280mm cuff length)	4690	Small	
	4691	Medium	
	4692	Large	
	4693	Extra Large	
Sleeve Extenders (Neoprene)	BNS-18-E	One	12 pair per bag/8 bags per carton
Sleeve Extenders (Vinyl)	VSC-HS-L	One	12 pair per bag/18 bags per carton
Sleeve Extenders (Polyethylene)	59-201	One	2000 pieces per carton



Ansell Limited
 Level 3, 678 Victoria Street
 Richmond, VIC 3121 Australia
 Telephone 61 3 9270 7270
 Facsimile 61 3 9270 7300
 protection@ap.ansell.com
 www.professional.ansell.com.au

Ansell and Micro-Touch® NitraTex® are trademarks owned by Ansell Limited, or one of its affiliates. © 2007 All Rights Reserved.

Ansell