



MICROTOUCH[®] DENTA-GLOVE[®]

BLUE NITRILE

PERFORMANCE SERIES



STAY COMFORTABLE AND RELAXED WITHOUT LOSING FUNCTIONALITY

Let your gloves work with your hands, not against them. Our MICROTOUCH DENTA-GLOVE Blue Nitrile gloves feature our innovative ERGOFORM™ ergonomic design technology to help reduce muscle effort that can cause stress on joints, ligaments and tendons. Durable enough for rugged donning while providing great tactile sensitivity and grip for fine, detailed procedures. The perfect go-to choice for users looking to avoid natural rubber latex Type I allergies while not compromising comfort. You'll never get tired of this!

KEY FEATURES & BENEFITS

- Innovative film technology provides comfort and grip
- Robust strength helps avoid rips and tears while donning
- ERGOFORM™ ergonomic design technology reduces hand fatigue during long procedures

RECOMMENDED FOR

- Enhanced tactile sensitivity and dexterity, especially for delicate procedures
- Standard, moderate risk examination procedures
- Extra comfort during long, sustained procedures
- Protection from latex Type I allergy in HCW's or patients
- Oral examinations

FEATURED TECHNOLOGIES



MICROTOUCH[®] premium medical examination gloves are scientifically designed to complement every need of medical professionals. Sophisticated design features include Ansell ergonomic technologies that ensure comfort during extended procedures, wet grip innovations that prevent harmful mishandling, and micro-thin films for maximum sensitivity. Caregivers who are fully prepared and comfortably protected by versatile MICROTOUCH styles have the peace of mind to focus on what matters most: the health and protection of themselves and their patients.



MICROTOUCH[®] DENTA-GLOVE[®]

BLUE NITRILE

PERFORMANCE SERIES



GENERAL INFORMATION

Material	Nitrile
Color	Blue
Certifications	In Compliance with ASTM D6319, D412, D5151, D6124, EN 455-1, EN 455-2, EN 455-3, EN 455-4, and ANSI/ASQ Z1.4 Manufacturing facility registered ISO 9001, ISO 13485
Powder content	Powder-free
Country of Origin	Malaysia

PHYSICAL PROPERTIES

	Before Aging	After Aging
Freedom from holes (inspection level 1)	1.5 AQL	N/A
Tensile, MPa (psi)	≥ 17	≥ 17
Elongation (%)	≥ 700	≥ 650
Force at Break (N)	≥ 6N	≥ 6N
Thickness (Single)	Typical Average Values	
Fingertip, mm (mils)	0.11 (4.3)	
Palm mm (mils)	0.08 (3.1)	
Cuff mm (mils)	0.07 (2.8)	

PACKAGING & STORAGE

Packaging	100 gloves/box, XL 90 gloves/box - 10 boxes/carton
Dispenser Dimensions	240mm x 130mm x 63mm

DIMENSIONAL SPECIFICATIONS

Minimum Length mm (inches)

X-Small	240 (9.4)
Small	240 (9.4)
Medium	240 (9.4)
Large	240 (9.4)
X-Large	240 (9.4)

REORDER (Size & Code)

X-Small	3130050XS
Small	31300500S
Medium	31300500M
Large	31300500L
X-Large	3130050XL

➤ **Contact Ansell Healthcare or your preferred distributor for a free sample.**

North America
US Tel: 800 952 9916
CA Tel: 1-844-494-7854
Email: insidesalesus@ansell.com

Europe, Middle East & Africa
Tel: +32 (0) 2 528 74 00
Email: info@ansell.eu

Asia Pacific
Tel: +603 8310 6688

Australia & New Zealand
Tel: 1800 337 041
Email: protectionau@ansell.com