



# MICROTOUCH<sup>®</sup> DENTA-GLOVE<sup>®</sup>

GREEN NEOPRENE

PERFORMANCE SERIES



## EVERYTHING YOU NEED IN ONE EXCEPTIONAL GLOVE

A bold claim that won't disappoint! We designed MICROTOUCH DENTA-GLOVE Green Neoprene to deliver the ultimate fit, comfort, tactile sensitivity, stretchability and grip. Make the switch from your current latex or nitrile and get the best performance without the limitations of each material. Featuring our ERGOFORM™ ergonomic design technology, this glove helps reduce muscle effort that can cause stress on joints, ligaments and tendons. Make the switch – your hands won't be disappointed!

### KEY FEATURES & BENEFITS

- Unique formulation for extra sensitivity with flexibility and comfort
- Strong and durable to avoid snagging and tearing through long procedures and donning
- Low to no odor
- ERGOFORM™ ergonomic design technology reduces hand fatigue during long procedures

### RECOMMENDED FOR

- Enhanced grip in wet environments
- Enhanced tactile sensitivity and dexterity, especially for delicate procedures
- Standard, moderate risk examination procedures
- Extra comfort during long, sustained procedures
- Protection from latex Type I allergy
- Oral examinations

### FEATURED TECHNOLOGIES



MICROTOUCH<sup>®</sup> premium medical examination gloves are scientifically designed to complement every need of medical professionals. Sophisticated design features include Ansell ergonomic technologies that ensure comfort during extended procedures, wet grip innovations that prevent harmful mishandling, and micro-thin films for maximum sensitivity. Caregivers who are fully prepared and comfortably protected by versatile MICROTOUCH styles have the peace of mind to focus on what matters most: the health and protection of themselves and their patients.



# MICROTOUCH<sup>®</sup> DENTA-GLOVE<sup>®</sup>

GREEN NEOPRENE

PERFORMANCE SERIES



## GENERAL INFORMATION

<b>Material</b>	Neoprene
<b>Color</b>	Green
<b>Certifications</b>	In Compliance with ASTM D6977, D412, D5151, D6124, EN 455-1, EN 455-2, EN 455-3, EN 455-4, and ANSI/ASQ Z1.4 Manufacturing facility registered ISO 9001, ISO 13485
<b>Powder content</b>	Powder-free
<b>Country of Origin</b>	Malaysia

## PHYSICAL PROPERTIES

	Before Aging	After Aging
<b>Freedom from holes</b> (inspection level 1)	1.5 AQL	N/A
<b>Tensile, MPa (psi)</b>	≥ 17	≥ 17
<b>Elongation (%)</b>	≥ 700	≥ 650
<b>Force at Break (N)</b>	≥ 6N	≥ 6N
<b>Thickness (Single)</b>	Typical Average Values	
<b>Fingertip, mm (mils)</b>	0.11 (4.3)	
<b>Palm mm (mils)</b>	0.08 (3.1)	
<b>Cuff mm (mils)</b>	0.07 (2.8)	

## PACKAGING & STORAGE

<b>Packaging</b>	100 gloves/box, XL 90 gloves/box - 10 boxes/carton
<b>Dispenser Dimensions</b>	240mm x 130mm x 63mm

## DIMENSIONAL SPECIFICATIONS

### Minimum Length mm (inches)

<b>X-Small</b>	240 (9.4)
<b>Small</b>	240 (9.4)
<b>Medium</b>	240 (9.4)
<b>Large</b>	240 (9.4)
<b>X-Large</b>	240 (9.4)

## REORDER (Size & Code)

<b>X-Small</b>	3130060XS
<b>Small</b>	31300600S
<b>Medium</b>	31300600M
<b>Large</b>	31300600L
<b>X-Large</b>	3130060XL

➤ **Contact Ansell Healthcare or your preferred distributor for a free sample.**

**North America**  
US Tel: 800 952 9916  
CA Tel: 1-844-494-7854  
Email: insidesalesus@ansell.com

**Europe, Middle East & Africa**  
Tel: +32 (0) 2 528 74 00  
Email: info@ansell.eu

**Asia Pacific**  
Tel: +603 8310 6688

**Australia & New Zealand**  
Tel: 1800 337 041  
Email: protectionau@ansell.com