



MICROTOUCH® DENTA-GLOVE®

GREEN NEOPRENE

PERFORMANCE SERIES



EVERYTHING YOU NEED IN ONE EXCEPTIONAL GLOVE

A bold claim that won't disappoint! We designed MICROTOUCH DENTA-GLOVE Green Neoprene to deliver the ultimate fit, comfort, tactile sensitivity, stretchability and grip. Make the switch from your current latex or nitrile and get the best performance without the limitations of each material. Featuring our ERGOFORM™ ergonomic design technology, this glove helps reduce muscle effort that can cause stress on joints, ligaments and tendons. Make the switch – your hands won't be disappointed!

KEY FEATURES & BENEFITS

- Unique formulation for extra sensitivity with flexibility and comfort
- Strong and durable to avoid snagging and tearing through long procedures and donning
- Low to no odor
- ERGOFORM™ ergonomic design technology reduces hand fatigue during long procedures

RECOMMENDED FOR

- Enhanced grip in wet environments
- Enhanced tactile sensitivity and dexterity, especially for delicate procedures
- Standard, moderate risk examination procedures
- Extra comfort during long, sustained procedures
- Protection from latex Type I allergy
- Oral examinations

FEATURED TECHNOLOGIES



MICROTOUCH® premium medical examination gloves are scientifically designed to complement every need of medical professionals. Sophisticated design features include Ansell ergonomic technologies that ensure comfort during extended procedures, wet grip innovations that prevent harmful mishandling, and micro-thin films for maximum sensitivity. Caregivers who are fully prepared and comfortably protected by versatile MICROTOUCH styles have the peace of mind to focus on what matters most: the health and protection of themselves and their patients.



MICROTOUCH[®] DENTA-GLOVE[®]

GREEN NEOPRENE

PERFORMANCE SERIES



GENERAL INFORMATION

Material	Neoprene
Color	Green
Certifications	In Compliance with ASTM D6977, D412, D5151, D6124, EN 455-1, EN 455-2, EN 455-3, EN 455-4, and ANSI/ASQ Z1.4 Manufacturing facility registered ISO 9001, ISO 13485
Powder content	Powder-free
Country of Origin	Malaysia

PHYSICAL PROPERTIES

	Before Aging	After Aging
Freedom from holes (inspection level 1)	1.5 AQL	N/A
Tensile, MPa (psi)	≥ 17	≥ 17
Elongation (%)	≥ 700	≥ 650
Force at Break (N)	≥ 6N	≥ 6N
Thickness (Single)	Typical Average Values	
Fingertip, mm (mils)	0.11 (4.3)	
Palm mm (mils)	0.08 (3.1)	
Cuff mm (mils)	0.07 (2.8)	

PACKAGING & STORAGE

Packaging	100 gloves/box, XL 90 gloves/box - 10 boxes/carton
Dispenser Dimensions	240mm x 130mm x 63mm

DIMENSIONAL SPECIFICATIONS

Minimum Length mm (inches)

X-Small	240 (9.4)
Small	240 (9.4)
Medium	240 (9.4)
Large	240 (9.4)
X-Large	240 (9.4)

REORDER (Size & Code)

X-Small	3130060XS
Small	31300600S
Medium	31300600M
Large	31300600L
X-Large	3130060XL

➤ **Contact Ansell Healthcare or your preferred distributor for a free sample.**

North America
US Tel: 800 952 9916
CA Tel: 1-844-494-7854
Email: insidesalesus@ansell.com

Europe, Middle East & Africa
Tel: +32 (0) 2 528 74 00
Email: info@ansell.eu

Asia Pacific
Tel: +603 8310 6688

Australia & New Zealand
Tel: 1800 337 041
Email: protectionau@ansell.com